Polyclonal Rabbit Anti-Human Alpha-1-Fetoprotein
Code No./ Code/ Code-Nr. A 0008
 Edition/ Ausgabe 01.12.02

ENGLISH	
Intended use	For in vitro diagnostic use. Polyclonal Rabbit Anti-Human Alpha-1-Fetoprotein (AFP), code No. A 0008, is intended for the determination of human alpha-1-fetoprotein in immunocytochemistry (1, 2, 3), ELISA (4) and gel immunoprecipitation tech-niques. Positive results may aid in the identification of non-neoplastic and neoplastic liver disease, yolk sac tumours and mixed germ cell tumours. However, it is important that the antibody is not used alone, but together with a panel of antibodies. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
Introduction	AFP is a 70 kDa glycoprotein containing 590 amino acid residues. Cells of the embryonic yolk sac, fetal liver, and fetal intestinal tract synthesize this glycoprotein (5, 6). AFP is a very useful marker for the detection of numerous developmental defects and various pathologies, whether cancerous or not (7). Expression of AFP has been demonstrated in about 44% of hepatocellular carcinomas, and in gonadal and extragonadal germ cells tumours, including yolk sac tumours. AFP has not been detected in normal adult tissue by immuno-staining (5). However, traces of AFP are found in normal adult sera and in greater concentrations in maternal and foetal sera and amniotic fluid (8).
Reagent provided	Purified immunoglobulin fraction of rabbit antiserum in liquid form. In 0.1 mol/L NaCl, 15 mmol/L NaN ₃ . <u>Protein concentration g/L:</u> See label on vial. The titre variation between different lots of A 0008 is less than 10%. This is achieved by adjusting the titre of each individual lot to match the titre of an antibody reference preparation kept at -80 °C.
Immunogen	Alpha-1-fetoprotein isolated from human cord serum.
Specificity	The antibody reacts with human alpha-1-fetoprotein. Traces of contaminating antibodies have been removed by solid-phase absorption with human plasma proteins. The specificity of the antibody has been ascertained as follows: <u>Crossed immunoelectrophoresis:</u> Only the alpha-1-fetoprotein precipitation arch appears when using 12.5 µL A 0008 per square cm gel area against 2 µL of a human cord serum fraction. No precipitation arch is seen when the antibody is tested against 2 µL human plasma. Staining: Coomassie Brilliant Blue. <u>ELISA:</u> In indirect ELISA no significant reaction is observed with normal human plasma. As demonstrated by rocket immunoelectrophoresis, the antibody cross-reacts with the AFP-equivalent protein in cat and wildebeest (9).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. The product may be used in different techniques and in combination with different sample types and materials, therefore each individual laboratory should validate the test system applied.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, relevant controls should be run simultaneously with patient specimens. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	IMMUNOCYTOCHEMISTRY <u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with proteinase K or heat-induced epitope retrieval is not required. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.
Staining procedure	<u>Dilution:</u> Polyclonal Rabbit Anti-Human Alpha-1-Fetoprotein, code No. A 0008, may be used at a dilution range of 1:400-1:800 when applied on formalin-fixed, paraffin-embedded sections of human foetal liver, and using 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), code No. X 0936, diluted to the same protein concentration as the primary antibody. Unless the stability in the actual test system has been established, it is recommended to dilute the product immediately before use or dilute in DakoCytomation Antibody Diluent, code No. S 0809. <u>Visualization:</u> DAKO LSAB™+HRP kit, code No. K 0679, and DAKO EnVision™+HRP kits, code Nos. K 4008 and K 4010, are recommended. Follow the procedure enclosed with the selected visualization kit.
Performance characteristics	Cells labelled by the antibody display a cytoplasmic staining pattern (1). <u>Normal tissues:</u> In a 20-day old human conceptus, the antibody labels the endodermal cells of the inner surface of the yolk sac. In the amnion AFP is also demonstrated (10). <u>Abnormal tissues:</u> 189 orchidectomy specimens with germ cell tumours (95 pure seminomas and 94 non-seminomas) were studied with the antibody. All seminomas were AFP negative, while 66% of the non-seminomas were positively stained for AFP. 21% of pure embryonal carcinomas (EC), 25% of EC components in mixed tumours, all yolk sac tumour components, 20% of pure teratomas (T) and 47% of T components were AFP positive. Carcinoma-in-situ present in seminiferous tubules adjacent to the tumours was never positively stained for AFP (1). In stomach cancers with elevated serum AFP levels, 19/35 (54%) were labelled by the antibody (2).
(102267-001)	A 0008/EFG/SS/01.12.02 p. 1/4

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

References/ Références/ Literatur	<ol style="list-style-type: none">Jacobsen GK, Jacobsen M. Alpha-fetoprotein (AFP) and human chorionic gonadotropin (HCG) in testicular germ cell tumours. A prospective immunohistochemical study. Acta Path Microbiol Immunol Scand Sect A 1983;91:165-76. Kodama T, Kameya T, Hirota T, Shimosato Y, Ohkura H, Mukojima T, et al. Production of alpha-fetoprotein, normal serum proteins, and human chorionic gonadotropin in stomach cancer: Histologic and immunohistochemical analyses of 35 cases. Cancer 1981;48:1647-55. Perrone T, Steeper TA, Dehner LP. Alpha-fetoprotein localization in pure ovarian teratoma. An immunohistochemical study of 12 cases. Am J Clin Pathol 1987;88:713-17. MacDonald DJ, Kelly AM. The rapid quantitation of serum alpha-fetoprotein by two-site micro enzyme immunoassay. Clin Chim Acta 1978;87:367-72. Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology; London: Oxford University Press; 1999. p. 9-10. Mizejewski GJ. <i>α</i>-fetoprotein as a biologic response modifier: Relevance to domain and subdomain structure (minireview). Proc Soc Exp Biol Med 1997;215:333-62. Nunez EA. Biological role of alpha-fetoprotein in the endocrinologic field: Data and hypotheses. Tumor Biol 1994;15:63-72. Tietz NW. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders Company; 1995. p. 236-38. Hau J, Nilsson M, Skovgaard-Jensen H-J, de Souza A, Eriksen E, Wandall LT. Analysis of animal serum proteins using antisera against human analogous proteins. Scand J Lab Anim Sci 1990;17:3-7. Jacobsen GK, Jacobsen M, Henriksen OB. An immunohistochemical study of a series of plasma proteins in the early human conceptus. Oncodev Biol Med 1981;2:399-410.		
Explanation of symbols/ Légende des symboles/ Erläuterung der Symbole			
REF	Catalogue number Référence du catalogue Bestellnummer	 <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p>	 <p>Manufacturer Fabricant Hersteller</p>
IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	LOT	Batch code Code du Lot Chargenbezeichnung
	Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten		Use by Utiliser jusque Verwendbar bis

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