**Artisan™  
Gram Stain Kit**

---

**ENGLISH  
Code AR175****Intended use**

For In Vitro Diagnostic Use.

The Gram Stain Kit is intended for laboratory use to identify, by light microscopy, Gram-positive and Gram-negative microorganisms in tissue samples.<sup>1-3</sup> Routinely processed samples (paraffin-embedded) may be used. The preferred fixative is neutral buffered formalin. The clinical interpretation of any staining, or its absence, should be complemented by morphological studies and proper controls, and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. All reagents in the Gram Stain Kit have been optimally prepared for use on Dako's Artisan Staining System and require no mixing or diluting.

**Summary and explanation**

The Gram Stain Kit is used to identify two distinct groups of microorganisms in tissue sections. Those that retain primary dye (Crystal Violet) are called Gram-positive. Those that lose the primary dye during decolorization are called Gram-negative.<sup>3</sup> Mechanisms of Gram-positive organisms retaining the primary stain and Gram-negative organisms losing the primary stain are the result of the chemistry and structure of the organism's cell walls.<sup>1</sup> This procedure is based on the work of Twort,<sup>1</sup> and utilizes a stable Crystal Violet-Ammonium Oxalate solution that aids in differentiation.

**Principle of procedure**

Bacteria are partially visible in Hematoxylin and Eosin preparations using high-dry and oil immersion light microscopy. However, greater sensitivity to the presence and specificity in their identification is achieved with the application of this special stain. The Gram Stain Kit is a modification of the original Twort's procedure substituting alcohol for acetone in the differentiating step.<sup>1</sup>

The Gram Stain Kit is optimized for use on the Artisan Staining System. Wash steps follow all of the staining steps. Following staining, slides are removed from the instrument, cleared with xylene or xylene substitute and mounted with appropriate mounting media. Results are interpreted using a standard light microscope.

**Controls:** Tissue containing both Gram positive and Gram negative bacteria.

**Reagents provided**

The Gram Stain Kit is composed of four ready-to-use dispenser packs, containing:

**Crystal Violet****CRYSTAL VIOLET**

2.3% Crystal Violet, 0.1% Ammonium Oxalate, 19% Ethyl Alcohol and 1% Methanol in deionized water containing an antimicrobial agent

**Lugol's Iodine****LUGOL'S IODINE**

1% Iodine and 2% Potassium Iodide in deionized water

**Twort's Fast Green****TWORT'S FAST GREEN**

0.2% Fast Green, 95% Alcohol

**Twort's Neutral Red****TWORT'S NEUTRAL RED**

0.2% Neutral Red, 95% Ethyl Alcohol, and 5% Methanol in deionized water

**The Crystal Violet reagent pack should be agitated before priming.**

Dako provides these reagents at optimal formulation for use according to the instructions provided for staining on prepared tissue sections. Any deviation from the recommended test procedures may invalidate declared expected results.

**Materials required, but not supplied**

1. Microscope slides
2. Artisan Staining System
3. Wash Solution (code AR102)

**Precautions**

1. For professional users.
2. Minimize microbial contamination of reagents or increase in non-specific staining may occur.
3. Reagent packs may contain residual solution after completing the prescribed number of tests.

4. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any changes made to the factory-released procedure.
5. Store reagent packs in an upright position when not in use on the instrument.
6. Keep reagent packs out of direct sunlight and away from heat-generating sources.
7. As a general rule, persons under 18 years of age are not allowed to work with this product. Users must be carefully instructed in the proper work procedures, the dangerous properties of the product and the necessary safety instructions. (Per European Union Directive 94/33/EC). Please refer to Material Safety Data Sheet (MSDS) for additional information.
8. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
9. Unused solution should be disposed of according to local, State and Federal regulations.
10. Safety Data Sheet available for professional users on request.

#### Risk and Safety statements

**Crystal Violet:** 3-<10% Methanol / Hazard Symbol: Harmful

R20/21/22 Harmful by inhalation, in contact with skin and if swallowed.

R68/20/21/22 Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed.

S35 This material and its container must be disposed of in a safe way.

S36/37 Wear suitable protective clothing and gloves.

**Lugol's Iodine:** 1-5% Potassium Iodide / Hazard Symbol: Irritant

R43 May cause sensitization by skin contact.

S60 This material and/or its container must be disposed of as hazardous waste.

**Twort's Fast Green:** 60-100% Ethanol, 3-<10% Methanol / Hazard Symbol: Highly Flammable, Harmful

R11 Highly flammable.

R20/21/22 Harmful by inhalation, in contact with skin and if swallowed.

R68/20/21/22 Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed.

S16 Keep away from sources of ignition – No smoking.

S35 This material and its container must be disposed of in a safe way.

S36/37 Wear suitable protective clothing and gloves.

**Twort's Neutral Red:** 60-100% Ethanol, 3-<10% Methanol / Hazard Symbol: Highly Flammable, Harmful

R10 Flammable.

S35 This material and its container must be disposed of in a safe way.

#### Storage

Each component should be stored at the temperature indicated on the label. Do not use after expiration date. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, appropriate tissue must be run as controls for such verification. If unexpected staining is observed which cannot be explained by variations in laboratory procedures or improper storage, contact Dako Technical Support.

#### Instrument procedure/ staining procedure

##### *Preliminary Comments*

Please refer to the *Artisan User Guide* for detailed instructions.

##### *Priming Reagent Packs*

Before beginning staining, reagent packs should be primed and checked for leakage. Check the prime of each reagent pack when a new kit is used and before the first staining run of the day.

1. Over a paper towel, if present, remove the red or black shipping cap from the dispense assembly of the reagent pack by gently pulling down on the cap.
2. To prime a reagent pack when it is first received, hold the reagent pack with the dispenser assembly pointed downward and lightly tap on the side of the reagent pack 2–3 times to release any air bubbles that may be trapped in the dispenser assembly.
3. Fully compress and release the plunger 3–4 times in rapid succession, taking care to absorb any expelled fluid into a paper towel, absorbent gauze, or similar material, until a steady stream is obtained.
4. Avoid dispensing excess amounts of reagent during priming. Excessive priming will result in reagent loss, which may cause the reagent pack to run out of fluid before the stated number of tests is performed.
5. Before each use, reagent packs may be primed 1–2 times, checking for a steady stream.
6. Holding the reagent pack in one hand, slide the mounting guide on the back of the reagent pack into the slots on the reagent carousel until the pack is completely seated.

##### *Procedure*

Refer to the *Procedure Reports* section in the Artisan Host Software for the recommended detailed procedure.

After staining is complete, slides will be dry. Dip slides in xylene or xylene substitute and permanently mount.

**Note:** Fresh 100% Alcohol must be used in the bulk liquid container when using this kit. Do not dip stained slides in alcohols as this will further differentiate the stain.

#### Staining interpretation

Gram-positive organisms:..... **Blue**

Gram-negative organisms: ..... **Red**

Background:..... **Varying shades of blue/green**

## Performance characteristics

### Reproducibility

Intra-run reproducibility of 100% was determined by staining 12 slides containing the same control tissue on one instrument.  
Inter-run reproducibility of 100% was determined by staining 12 slides containing the same control tissue on three different instruments.

## Troubleshooting

Refer to the Maintenance and Troubleshooting section of your *Artisan User Guide* or contact our Technical Service Group.

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#### Bibliographie

#### Literaturhinweise

1. Bancroft JD, Stevens A. Theory and Practice of Histological Techniques, 4th edition, Churchill Livingstone Inc., New York, New York, 1996
2. Luna LG. Manual of Histologic Staining Methods of the AFIP, 3rd edition, McGraw-Hill Book Company, New York, New York, 1968
3. Prophet EB, Mills B, Arrington JB, Sobin LH. AFIP: Laboratory Methods in Histotechnology, 1st edition, American Registry of Pathology, Washington, DC, 1992

Highly flammable Extrêmement inflammable Leichtentzündlich PT0062 Rev A		Harmful Nocif Gesundheitsschädlich PT0061 Rev A		Irritant Irritant Reizend PT0060 Rev A	
Catalogue number Référence du catalogue Bestellnummer PT0020/ Rev C	Temperature limitation Limites de température Zulässiger Temperaturbereich	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum			
Manufacturer Fabricant Hersteller	Batch code Code du lot Chargenbezeichnung	Contains sufficient for <N> tests Contenu suffisant pour <N> tests Inhalt ausreichend für „N“ Ansätze	Dako North America, Inc. 6392 Via Real Carpinteria, California 93013 USA		
Use by Utiliser jusque Verwendbar bis	Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	Tel +45 4485 6655 Fax 805 566 6688 Technical Support 800 424 0021 Customer Service 800 235 5763 www.dako.com		
			Dako Denmark A/S Produktionsvej 42 DK-2600 Glostrup Denmark Tel +45 4485 9500 Fax +45 4485 9595		

Edition 08/07