**Monoclonal Mouse
Anti-Human
N-Cadherin
Clone 6G11**

ENGLISH
Code M3613

Intended use

For In Vitro Diagnostic Use.

This antibody is intended for laboratory use to identify qualitatively by light microscopy N-cadherin expressing cells in normal and neoplastic tissues using immunohistochemical (IHC) test methods. Positive results aid in the classification of mesotheliomas.²⁻⁴ Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Synonym

Nerve cadherin, Neural-type cadherin (NCAD)

Summary and explanation

N-cadherin is a 140 kD protein that belongs to a family of transmembrane molecules that mediate calcium-dependent intercellular adhesion. Cadherins are involved in controlling morphogenetic movements during development and regulate cell surface adhesion through homotypic adhesion with the same cadherin species. N-cadherin's function is dependent on its association with the actin-cytoskeleton, which is mediated through interactions between the C-terminal region of N-cadherin and the cytoplasmic catenin proteins. The stability of this association is regulated by phosphorylation and dephosphorylation of β -catenin.^{5,6}

Expression of N-cadherin has been reported on a variety of normal tissues including neuronal, endothelial and muscle cells and a subpopulation of early hematopoietic progenitor cells.⁷ Among neoplastic tissues, N-cadherin expression has been demonstrated in malignant noncarcinomatous neoplasms including mesotheliomas,^{2-4,8} chordomas, synovial sarcomas, malignant melanomas, epithelioid sarcomas, epithelioid angiosarcomas and clear cell sarcomas.⁸ Serous and endometrioid tumors of the ovary have been demonstrated to be N-cadherin positive whereas mucinous tumors are negative.⁹ Other N-cadherin positive neoplasms include renal cell carcinomas and some variant breast tumors including medullary breast carcinomas and sarcomatoid metaplastic breast carcinomas.^{10,11}

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Anti-N-cadherin is provided as tissue culture supernatant (containing fetal bovine serum) dialyzed against 0.05 mol/L Tris-HCl, pH 7.2, and 0.015 mol/L sodium azide. Contains stabilizing protein.

Clone: 6G11¹ Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

M3613 may be used at a dilution 1:50 when performing IHC using the EnVision+™ detection system. These are guidelines only. Optimal antibody concentrations may vary depending on specimen and preparation method, and should be determined individually in each laboratory.

Immunogen

Synthetic peptide corresponding to human N-cadherin amino acids 878-905¹

Specificity

Clone 6G11 is specific to the peptide used for immunization and is unreactive with the KLH carrier protein. The sequence used for immunization is homologous among human, mouse, rat and chicken N-cadherin proteins. The sequence has no homology to other cadherin proteins such as E-cadherin, P-cadherin, Cadherin-6, Cadherin-12 or to the catenins (alpha, beta or gamma-catenin).¹

Materials required, but not supplied

Refer to the "General Instructions for Immunohistochemical Staining" and/or the Detection System "Instructions". Suggested diluent(s) for IHC procedures:

Antibody Diluent (code S0809)

The following negative control is recommended for IHC procedures:

Mouse IgG₁ (code X0931)

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin Sections

Anti-N-Cadherin can be used on formalin-fixed, paraffin-embedded tissue sections. Pretreatment of tissue with proteolytic enzymes is not recommended.

The deparaffinized tissue sections must be treated with heat prior to the IHC staining procedure. Heat-induced epitope retrieval (HIER) involves immersion of tissue sections in a pre-heated buffer solution and maintaining heat, either in a water bath or steamer (95–99 °C), or a pressure cooker (121 °C). For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended. Target Retrieval Solution (code S1700) or 10x Concentrate (code S1699) is recommended using a 40-minute heating protocol for HIER performed at 95–99 °C. After thermal treatment, allow the jar with buffer and slides to cool for 20 minutes at room temperature. Rinse well with buffer or deionized water. Optimal HIER incubation times at temperatures greater than 99 °C should be determined by each individual laboratory.

Cryostat Sections And Cell Smears

Anti-N-Cadherin can be used for labelling acetone-fixed cryostat sections or fixed cell smears.

Prior to IHC staining, tissues must be fixed and processed. Fixation prevents autolysis and putrefaction of excised tissues, preserves antigenicity, enhances the refractive index of tissue constituents and increases the resistance of cellular elements to tissue processing. Tissue processing includes dehydration, clearing of dehydrating agents, infiltration of embedding media, embedding and sectioning of tissues. The most common fixatives for IHC tissue preparations are discussed in the "General Instructions for Immunohistochemical Staining". These are guidelines only. Optimal procedures must be determined and verified by the user.

Staining procedure

Follow the recommended procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-N-cadherin is membrane and cytoplasmic.

Performance characteristics

Normal tissues¹²

Anti-N-cadherin strongly immunostains ganglion cells and nerve in many normal tissue types.

Tissue Type (# tested)	Positively Staining Tissue Elements
Adrenal (3)	3/3 membrane and cytoplasm of adrenal cortical and adrenal medulla cells
Brain, cerebellum (3)	3/3 diffuse staining of neuropil and membrane of Purkinje cells
Brain, cerebrum (3)	3/3 diffuse staining of neuropil
Stomach (3)	3/3 membrane and cytoplasm of chief, parietal and oxyphil cells
Endometrium (3)	3/3 membrane and cytoplasm of endometrial glands and cytoplasm of endometrial stromal cells
Kidney (3)	3/3 membrane and cytoplasm of renal tubule epithelium
Ovary (3)	1/3 membrane and cytoplasm of surface epithelial cells and inclusion cysts epithelium
Pancreas (3)	3/3 membrane and cytoplasm of acinar cells
Pituitary (3)	3/3 diffuse staining of cell processes in neurohypophysis and membrane and cytoplasm of pituicytes in adenohypophysis
Heart (3)	3/3 cell junction and cytoplasm of 100% of myocytes
Liver (3)	3/3 membrane and cytoplasm of hepatocytes and membrane of bile duct epithelium
Spleen (3)	3/3 membrane and cytoplasm of splenic sinuses
Testis (3)	3/3 membrane and cytoplasm of 100% of sertoli cells, spermatogonia and seminiferous tubules
Tonsil (3)	3/3 cytoplasm germinal center cells
Parathyroid (3)	3/3 membrane and cytoplasm of chief cells
Skeletal Muscle (3)	0/3
Colon (3)	0/3
Esophagus (3)	0/3
Lung (3)	0/3
Skin (3)	0/3
Bone Marrow (3)	0/3
Breast (3)	1/3 membrane of ductal epithelial cells (1 duct)
Cervix (3)	0/3
Prostate (3)	0/3
Thymus (3)	0/3
Mesothelial Cells (3)	0/3
Nerve, peripheral (3)	0/3

Salivary Gland (3)	0/3
Small Intestine (3)	0/3
Thyroid (3)	0/3

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Bibliographie

Literatur

- Report on File
- Abutally AS, Addis BJ, Roche WR. Immunohistochemistry in the distinction between malignant mesothelioma and pulmonary adenocarcinoma: a critical evaluation of new antibodies. J Clin Pathol. 2002; 55(9):662
- Thirkettle I, Harvey P, Hasleton PS, Ball RY, Warn RM. Immunoreactivity for cadherins, HGF/SF, met, and erbB-2 in pleural malignant mesotheliomas. Histopathology. 2000; 36(6):522
- Han AC, Peralta Soler A, Knudsen KA, Wheelock MJ, Johnson KR, Salazar H. Differential expression of N-cadherin in pleural mesotheliomas and E-cadherin in lung adenocarcinomas in formalin-fixed, paraffin-embedded tissues. Hum Pathol. 1997; 28(6):641
- Takeichi M. The cadherins: cell-cell adhesion molecules controlling animal morphogenesis. Development. 1988; 102(4):639
- Balsamo J, Arregui C, Leung T, Lilien J. The nonreceptor protein tyrosine phosphatase PTP1B binds to the cytoplasmic domain of N-cadherin and regulates the cadherin-actin linkage. J Cell Biol. 1998; 143(2):523
- Puch S, Armeanu S, Kibler C, Johnson KR, Muller CA, Wheelock MJ, Klein G. N-cadherin is developmentally regulated and functionally involved in early hematopoietic cell differentiation. J Cell Sci. 2001; 114(8):1567
- Laskin WB, Miettinen M. Epithelial-type and neural-type cadherin expression in malignant noncarcinomatous neoplasms with epithelioid features that involve the soft tissues. Arch Pathol Lab Med. 2002; 126(4):425
- Peralta Soler A, Knudsen KA, Tecson-Miguel A, McBrearty FX, Han AC, Salazar H. Expression of E-cadherin and N-cadherin in surface epithelial-stromal tumors of the ovary distinguishes mucinous from serous and endometrioid tumors. Hum Pathol 1997; 28(6):734
- Markovic-Lipkovski J, Brasanac D, Muller GA, Muller CA. Cadherins and integrins in renal cell carcinoma: an immunohistochemical study. Tumori. 2001;87(3):173
- Han AC, Peralta Soler A, Knudsen KA, Salazar H. Distinct cadherin profiles in special variant carcinomas and other tumors of the breast. Hum Pathol. 1999; 30(9):1035
- IHC003 Report On File

REF Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten
 Manufacturer Fabricant Hersteller	LOT Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis
EC REP Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	IVD In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA

EC REP

Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark

Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

Tel +45 4485 9500
Fax +45 4485 9595

www.dako.com

PT0039/Rev C